



### USER-FRIENDLY DIAGNOSTICS

- ◆ Icon driven touch screen is intuitive and is easy to navigate.
- ◆ Assay monitoring enables the user to check sample progress and assay performance
- ◆ The system can be interfaced to any LIS system
- ◆ Up to 300 assay results can be saved in memory



### ADVANCED SAMPLE CAROUSEL

- ◆ The AIA-360 sample carousel holds up to 25 patient samples and 25 test cups at one time
- ◆ Up to four assays per specimen can be loaded
- ◆ Additional samples and tests can be added at any time without interrupting



### FAST AND ACCURATE

- ◆ First result in ~20 minutes
- ◆ Bar-code sample ID ensures accurate patient identification without manual worklist entry
- ◆ On-board printer delivers patient results in real-time for fast and easy review
- ◆ Primary tube sampling



### CONVENIENT REAGENT ACCESS

- ◆ The common reagents are easily accessible from the removable reagent tray
- ◆ Reagent level sensors signal the user when reagents are running low or when the waste needs to be emptied



KOMPLETNO MEDICINSKO SNABDEVANJE  
Laze Lazarevića 25, 21000 Novi Sad, Srbija  
Tel / Fax: +381 / (021) 466-800; 466-809; 468-185



TOSOH Bioscience, Inc.  
6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080  
Tel: (800) 248-6764 • Fax: (650) 615-0415  
www.tosoh.com

### AIA-360 Result Printout

```

OP:SHANTI      05/11/04 08:59
SampleID:100010
Samp. Pos:1
Analyte :#FER  Lot:11 :1
Result  :79.2 ng/ml
Flag    :
Analyte :#bHCG Lot:33 :2
Result  :289.6 mIU/ml
Flag    :
    
```